



Pole mounting hardware pH, 10cm bracket, SS pole 2m

Nr. produs: LZY714.99.23920
RON Preț (fără TVA): Contact
Disponibil

Stainless steel pole mounting hardware for 1" NPT pH sensor including 10 cm SS bracket and 2 m SS pole. Suitable for all 1" NPT pH sensors like 1200S sc (LXV426.99.10001) and pHD ORP (LXV427.99.10001)

Specificații

Diametru:	40 mm
Greutate:	Lifting part only: 4.9 kg
Instrument:	Suitable for all 1" NPT pH sensors
Lungime:	Pole: 2000 mm
Material:	Wetted parts: pole SS 1.4301; sensor adapter PVC; screw set SS.