



Beaker, polypropylene, low form, 600 mL

Nr. produs: 108052
RON Preț (fără TVA): Contact
Expediere în 2 săptămâni

Griffin low form. For general laboratory use. Non-toxic, noncontaminating, and biologically inert. Autoclavable. Combines no drip pouring, durability, and translucency. Approximate volume scales molded in for small sizes. Sizes 250 mL and larger feature silk-screened graduations and double spout for easier pouring. Nalgene.

Specificații

Capacitate:	600 mL
Diametru exterior :	110.7 mm with spout
Instrument:	All except DR/4000, PC II
Material:	Polypropylene, Low Form
Numar metoda:	8006
Unitate:	1 each