



## Beaker, polypropylene, low form, 100 mL, 12/pk

**Nr. produs:** 108072  
**RON Preț (fără TVA):** Contact  
**Expediere în 2 săptămâni**

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes  $\geq 250$  mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

---

### Specificații

Capacitate:	100 mL
Diametru exterior :	63.0 mm with spout
Material:	Polypropylene, Low Form
Material corp:	Polypropylene
Unitate:	12 /pk