



Be Right™



Funnel, analytical, 65 mL approx. volume

Nr. produs: 108367

RON Preț (fără TVA): Contact

Expediere în 2 săptămâni

Standard long stem funnel meets the needs of a wide variety of applications. Mounded ribbing prevents "air lock" commonly occurring during funneling and filtrations. A steep, 58-degree cone slope makes for fast and easy flow of either liquids or small solids. Polypropylene construction is resistant to chemicals and makes this funnel fully autoclavable for sterilization.

Specificații

Cantitate: 1 each