



Reagent powder pillows, ChromaVer® 3, chromium, , 5 or 10 mL, pk/1000

Nr. produs:

1271028

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



For determination of hexavalent chromium by the Diphenylcarbohydrazide method. USEPA accepted for reporting for wastewater analyses using Hach Method 8023, Diphenylcarbohydrazide Method. Photometric Sample Size: 10 mL. Photometric Range: up to 0.70 mg/L. Kit Sample Size: 5 mL. Kit Range: Up to 1.50 Pack of 1000 powder pillows.

Specificații

Domeniu de măsurare:	0.010 - 0.700 mg/L Cr ⁶⁺
EPA compliant:	Da
Instrument:	All except DR/3000, DR/2000, DR 900
Numar metoda:	8023
Numar teste:	1000
Nume metoda:	1,5-Diphenylcarbohydrazide
Parametru:	Chromium
Platforma:	Pernu#e cu pudră