



Nitrate Standard Solution, 1000 mg/L, 500 mL

Nr. produs:

1279249

RON Preț (fără TVA):

Contact

Disponibil

Produse periculoase



Ion Selective Electrode (ISE) standard solutions for laboratory analyses

ISE Standard Solution. Nitrate, 1000 mg/L ± 10 mg/L $\text{NO}_3\text{-N}$ (NIST). 500 mL bottle

Traceable to NIST standard reference materials

Always performing measurement results you can trust is essential. The quality of the standard solutions that are used are key.

Design for use with Hach meters portfolio and other meters.

The basis for accurate ISE measurements is good care and calibration of the Ion Selective Electrode.

Reliable packaging

Supplied in an airtight packaging ensuring guaranteed shelf-life.

Specificații

Cantitate:	500 mL
Concentrație:	4430 ± 44 mg/L NO_3
Condiții de depozitare:	10 - 25 °C (protect from light)
Descriere:	Nitrogen-Nitrate Standard Solution as $\text{NO}_3\text{-N}$ (NIST)
Note speciale:	Observed shelf life after opening, following recommended storage and use: 3 months
Package settings:	Bottle
Parametru:	Nitrate
Platforma:	Soluție