



Color standard solution, 500 Pt Co units, pk/16 - 10-mL Voluette® ampules

Nr. produs: 141410
RON Preț (fără TVA): Contact
Disponibil

Produse periculoase



For accuracy checks in Color determinations. 500 platinum-cobalt units. pk/16 - 10 mL Voluette® Ampules.

Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

Reliable results

Sealed ampules prevent accidental contamination or evaporation of standard solutions.

Specificații

Cantitate:	10 mL
Concentrație:	500 units
Descriere:	Color Standards
Platforma:	Solu#ie