



Digital titrator cartridge, hydrochloric acid , 8.00 N

Nr. produs:

1439001

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



For base concentration determination by Digital Titrator titration. Digits x 0.01 = meq

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specificații

Domeniu de măsurare:	10 - 4000 meq/L Base
Numar metoda:	8200 Acid)
Numar teste:	40 - 100
Nume metoda:	Phenolphthalein End Point
Parametru:	Acid-Base
Platforma:	Titrator digital
Titrant:	8.00 N HCl