



CDTA Digital titrator cartridge, 0.0800 M

Nr. produs: 1440201

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Cyclohexanediamine tetraacetic acid standard solution. For Calcium and Total Hardness determination by Digital Titrator titration. 0.0800 M. Digits x 0.1 = mg/L CaCO₃ in 100 mL. Digits x 0.0001 = mmole.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specificații

Domeniu de măsurare:	10 - 160 mg/L CaCO ₃
Nume metoda:	CDTA Titration
Platforma:	Titratore digital

Accesorii necesare

- Mănuși de protecție mărimea L, albastre, nitril, fără pudră, 50 de bucăți (Item SM743L)