



Iron and manganese color disc test kit, model IR-20

Nr. produs:

146300

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



When customers complain of stained and dingy laundry, this kit lets you do quick, on-the-spot tests for iron and manganese to see where the problem lies.

Method / Chemistry: Color disc colorimetry. Test: Iron/ Manganese. Range (mg/L): 0-5/ 0-3. Smallest I

Simple, fast, color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Aplicație:	Water Quality
Cei mai mici pași de incrementare:	0.2
Domeniu de măsurare:	0 - 3.0
Domeniu de masurare :	0 - 4,0 mg/L Fe
Model:	IR-20
Parametru:	Iron and Manganese
Platforma:	Disc culoare