



Manganese test kit, model MN-5

Nr. produs: 146700
RON Preț (fără TVA): Contact
Nu este afișată data expirării

Produse periculoase



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. An economical solution, this kit uses the simple, accurate color disc method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specificații

Domeniu de măsurare:	0 - 3.0 mg/L Mn
Interval de măsurare standard:	0 - 3.0 mg/L Mn
Model:	MN-5
Model test kit:	MN-5
Numar teste:	100
Nume metoda:	Color disc/Cold periodate
Parameter:	Manganese
Parametru:	Manganese - As Mn
Platforma:	Disc culoare
Smallest increment:	0.1
Tipul carcasei:	D