



Laboratory Bottles, Polypropylene, Autoclavable, Wide Mouth, 1 L, 6/pk

Nr. produs: 1472463

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Polypropylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121 °C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Specificații

Body Material:	Plastic
Cantitate:	6 /pk
Capacitate:	1000 mL
Material:	PP
Package size UOM:	6pk