



Chromium color disc test kit, model CH-8

Nr. produs: 183400

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Produse periculoase



Widely used as a corrosion inhibitor in cooling towers and engine cooling systems, chromium is also used in surface finishing, leather tanning, and other industries. The most widely used form, hexavalent chromium, is very toxic and strict limits are placed on discharges. Utilizing the color disc method to determine chromium (Range: 0.1-1.5 mg/L), this kit contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

| | |
|----------------------|--|
| Domeniu de măsurare: | 0 - 1.3 Cr ⁶⁺ |
| Model: | CH-8 |
| Model test kit: | CH-8 |
| Numar teste: | 100 |
| Nume metoda: | Color disc/ Diphenylcarbazide |
| Parameter: | Chromium Hexavalent |
| Parametru: | Chromium, low range - As hexavalent Cr |
| Platforma: | Disc culoare |
| Smallest increment: | 0.1 |
| Tipul carcasei: | D |