



Cyanuric acid test kit, model CY-3

Nr. produs: 185102

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

A stabilizing agent for chlorine, cyanuric acid is used in swimming pool treatment. Cyanuric Acid Test Kit measures the turbidity formed in proportion to Cyanuric Acid concentration, after the reagent is added. Contains 100 tests.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Aplicatie:	Pool and Spa
Domeniu de măsurare:	20 - 100 mg/L
Model:	CY-3
Model test kit:	CY-3
Numar teste:	100
Nume metoda:	Turbidimetry
Parametru:	Cyanuric Acid
Smallest increment:	20
Tipul carcasei:	D