



Bottle, BOD, 60 mL capacity, marked 30 mL graduation line

Nr. produs:	190902
RON Preț (fără TVA):	Contact
Expediere în 2 săptămâni	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 mL sizes).

Specificații

Body Material:	Borosilicate Glass
Cantitate:	1 each
Capacitate:	60 mL
Material:	Borosilicate Glass
Package size UOM:	1pk