



Test kit tannin/lignin

Nr. produs: 193701
RON Preț (fără TVA): Contact
 Nu este afișată data expirării

Produse periculoase



A product of decomposed plant material, tannin occurs in natural waters and is used for boiler and cooling tower water treatment, dyeing, ink manufacture, tanning, and paper sizing. Lignin, also a plant product, is found in paper pulp manufacturing waste. Tannin and lignin are analyzed together and reported as mg/L tannic acid. This kit utilizes the simple, low-cost color disc method - the most accurate visual test kit method.

Tannin and lignin are analyzed together and reported as mg/L tannic acid.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specificații

Domeniu de măsurare:	0 - 15 mg/L as Tannic Acid
Domeniu de masurare :	0 - 150 mg/L
Interval de măsurare standard:	0 - 15 mg/L as Tannic Acid; 0 - 150 mg/L
Model:	Tannin/Lignin
Model test kit:	TA-3
Numar teste:	100
Nume metoda:	Color disc/Tyrosine
Parameter:	Tannin/Lignin
Parametru:	Tannin/Lignin - As tannic acid
Platforma:	Disc culoare
Smallest increment:	0.5
	5
Tipul carcasei:	R