



Test kit, hardness (total), model HA-DT

Nr. produs:

2063600

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Utilizing the Hach Digital Titrator, the Model HA-DT kit is used to determine hardness concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-4,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specificații

Cei mai mici pași de incrementare:	0.1-10
Domeniu de măsurare:	10 - 4,000 CaCO ₃
Model test kit:	HA-DT
Numar teste:	100
Nume metoda:	Digital Titrator/EDTA
Parametru:	Hardness, total - As CaCO ₃
Platforma:	Titrator digital
Tipul carcasei:	F