



Tubing, nalgene 8000, vacuum, 9.5 mm (I.D.)

Nr. produs: 2074146

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Similar to Nalgene 8000 regular tubing, but with heavier walls. Imprinted in red for instant identification. Saves time and avoids errors which could invalidate experiments or injure laboratory personnel.

Specificații

Diametru exterior :	7/8 in
Formula (Text):	R-3603
Grosime conducte de perete:	1/4 in
I. D.:	3/8 in
Lungime:	50 ft
Material corp:	Tygon