



Bottle, Sampling, Low Density, Polyethylene, 250 mL, 12/pk

Nr. produs:	2087076
RON Preț (fără TVA):	Contact
Expediere în 2 săptămâni	

Low-density polyethylene (LDPE) bottles feature wide-mouth opening for easy filling and leak-proof polypropylene closure. Great choice for grab samples, shipping liquids, soil sampling, and storing dry chemicals.

Specificații

Body Material:	Low Density Plastic
Cantitate:	12 /pk
Capacitate:	250 mL
Material:	LDPE
Package size UOM:	12pk