



Cylinder, Graduated, Mixing, Glass, 10 mL ± 0.1 mL, 0.2 mL divisions, glass stopper #13

Nr. produs: 2088638

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

This glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes glass penny head stopper. Graduations are in durable white enamel bake.

Specificații

Cantitate:	1 each
Capacitate:	10 mL
Certificări:	Class B
Material corp:	Glass
Metoda de măsurare:	To Contain
Tip:	Simplu Metric, Mixing