



Be Right™



Cylinder, graduated, mixing, Glass, 50 mL ± 0.4 mL, 1.0 mL divisions, glass stopper #16

Nr. produs: 2088641

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

This tall form glass graduated cylinder is calibrated “to contain” with a molded in base. Each includes glass penny head stopper. Graduations are in durable white enamel bake.

Specificații

Cantitate:	1 each
Capacitate:	50 mL
Certificări:	Class B
Material corp:	Glass
Metoda de măsurare:	To Contain
Tip:	Simplu Metric, Mixing