



**Be Right™**



## Rack, drying

**Nr. produs:** 2089000

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Epoxy-coated, steel-wire drying rack for test tubes, vials, flasks, beakers, and other glassware. Three rows of hair-pin looped supports to hold large-diameter necks. Fourteen loops total. Top row has 12 straight-pin for smaller tubes and glassware. Polystyrene tray fits into bottom rails to catch draining water. Dimensions: 43 x 17 x 48 cm (17 x 6.7 x 19"). Free-standing or wall mount option.

---

### Specificații

Adancime:	19 Inches
Cantitate:	1 each
Lungime:	430 mm