



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 500 mL

Nr. produs: 2089649

ARTICOL ÎNVECHIT

Acest articol nu mai este disponibil.

Produse înlocuitoare sugerate:

- 2089679

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Specificații

Capacitate:	500 mL
Material:	Glass Screw Cap
Unitate:	1 each