



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 6/pk, 250 mL

Nr. produs:

2089676

ARTICOL ÎNVECHIT

Acest articol nu mai este disponibil.

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Specificații

Capacitate:	250 mL
Material:	Glass Screw Cap
Unitate:	6/pk