



Flask, Erlenmeyer, Polycarbonate Capacity 125 mL

Nr. produs: 2089743

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Graduated (Molded-in graduations) flask features: clarity, toughness and strength. Autoclavable; however, repeated autoclaving may result in loss of mechanical strength. #5 Neoprene stoppers (P/N 14808-05; pk/12) available separately. 125 mL.

Specificații

Capacitate:	125 mL
Material:	Polycarbonate
Unitate:	1 each