



Be Right™



Flask, Erlenmeyer, polymethylpentene capacity 125 mL, 6/pk

Nr. produs: 2089873

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

Specificații

Capacitate:	125 mL
Material:	Polymethylpentene
Tip:	Erlenmeyer
Unitate:	6/pk