



Pipet rack

Nr. produs:

2094700

RON Preț (fără TVA):

Contact

Nu este afișată data expedierii

Polypropylene, fifty-place rack for safe, orderly storage of pipets. Helps create an efficient storage system or a well organized liquid handling station. Holes in top and middle decks are 16 mm diameter (0.63"); bottom holes are tapered to 6 mm (0.24") to prevent contact with supporting surface. Autoclavable to 122°C. Overall dimension: 21 x 11 x 22 cm (8.3 x 4.3 x 8.7").

Specificații

Adâncime:	8.7 Inches
Cantitate:	1 each
Latime:	4.3 Inches
Lungime:	210 mm