



Cylinder, graduated, polymethylpentene, 250 mL \pm 1.4 mL, 2.0 mL divisions

Nr. produs: 217246
RON Preț (fără TVA): Contact
Nu este afișată data expirării

This cylinder features transparent, nontoxic, noncontaminating, and biologically inert material with excellent chemical resistance. Meets ASTM E1272 Class B accuracy requirements, ISO Standard 6706, and suitable for food and beverage use (PMP resin meets requirements of CFR 21, Part 177.1520). Calibrated "To Contain/To Deliver" at 20°C (68°F). Shatterproof and thermally resistant to 135°C (275°F). Absence of meniscus for easy-to-read levels on front-to-back raised graduations.

Specificații

Cantitate:	1 each
Capacitate:	250 mL
Certificări:	Class B
Material corp:	Polymethylpentene
Metoda de măsurare:	To Deliver/To Contain
Tip:	Simplu Metric