



Test kit, formaldehyde, model FM-1

Nr. produs:

2183100

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Formaldehyde solutions are used for maintenance of reverse osmosis systems and as disinfectants. The Model FM-1 measures higher concentrations by drop count titration. Contains 100 tests.

Utilises drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Domeniu de măsurare:	0,05 - 1% & 0,5 - 10%
Model test kit:	FM-1
Numar teste:	100
Nume metoda:	Drop count titration/
	Thymolphthalein
Parametru:	Formaldehyde - As CH_2O
Platforma:	Drop Count Titration
Smallest increment:	0.05%
	0.5%
Tipul carcasei:	D