



Be Right™



Pipet tips, Tensette pipette 1970010, 1.0-10.0 mL, non-sterile, pk/250

Nr. produs: 2199725

RON Preț (fără TVA): Contact

Expediere în 2 săptămâni

Perform accurate, economical pipetting with Hach's TenSette pipets. Single-channel, air-displacement pipets come in two models: 0.1-1.0 mL volumes (0.1 mL increments) or 1.0-10.0 mL (1.0 mL increments). Volume selection with a simple twist of the turret.

Specificații

Cantitate: 250 /pk