



Test kit, chromium, model CH-14

Nr. produs: 222702
RON Preț (fără TVA): Contact
 Nu este afișată data expirării

Produse periculoase



Widely used as a corrosion inhibitor in cooling towers and engine cooling systems, chromium is also used in surface finishing, leather tanning, and other industries. The most widely used form, hexavalent chromium, is very toxic and strict limits are placed on discharges. Utilizing the drop count titration method to determine chromium (high range), this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Domeniu de măsurare:	5 - 100 mg/L & 50 - 1000 mg/L Cr ⁶⁺
Model test kit:	CH-14
Numar teste:	100
Nume metoda:	Drop count titration/ Thiosulfate
Parametru:	Chromium, high range - As Hexavalent Cr
Platforma:	Drop Count Titration
Smallest increment:	5 50
Tipul carcasei:	D