



Bag, sterile, WhirlPak, 207 mL, 100/pk

Nr. produs: 2233199

RON Preț (fără TVA): Contact

Expediere în 2 săptămâni

Sampling and Handling

Whirl-Pak Sterile Bag Vol: 207 mL (7oz), Dim (W x L): 9.5 cm x 18 cm, (3.75 x 7") 3.0 mil thick

Whirl-Pak Sterile Bag Vol: 207 mL (7oz), Dim (W x L): 9.5 cm x 18 cm, (3.75 x 7") 3.0 mil thick

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection. "Write-on" surface accepts waterproof ink. 24 oz bag perfect for incubating petri dishes in water bath applications.

Specificații

Cantitate:	100
Capacitate:	207 mL