



Reagent and apparatus set, water in oil

Nr. produs:

2237300

Produse periculoase

ARTICOL ÎNVECHIT

Acest articol nu mai este disponibil.



Model WO-1 Test Kit

Model WO-1 Test Kit

The Model WO-1 test kit provides critical data in oil quality determinations, oil line operation, and diesel engine maintenance. The kit has proven particularly useful on board ships. Reagent is packaged in unit-dose vials, and parts in contact with the sample are disposable, so clean-up is eliminated.

Measures low levels of water in oil

Provides critical data in oil quality deterinations, oil line operation, and diesel engine maintenance

Reagent is packaged in unit-dose vials, and parts in contact with the sample are disposable, so clean-up is eliminated

Specificații

Domeniu de măsurare:	0.05 - 1 %
Domeniu de masurare :	0 - 10 %
Model test kit:	WO-1
Numar teste:	25
Nume metoda:	Volumetric displacement/ Calcium hydride
Parametru:	Water-in-oil - As %
Smallest increment:	0.05% 0.20%
Tipul carcasei:	no case