



Test kit, nitrogen , model NI-8

Nr. produs: 224100
RON Preț (fără TVA): Contact
 Nu este afișată data expirării

Produse periculoase



Designed for testing ammonia in freshwater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

One of Hach's most popular test kits

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specificații

Domeniu de măsurare:	0 - 3.0 mg/L NH ₃ -N
Model test kit:	NI-8
Numar teste:	100
Nume metoda:	Color disc/Nessler reagent
Parametru:	Nitrogen, Ammonia, mid range - For freshwater
Platforma:	Disc culoare
Smallest increment:	0.1
Tipul carcasei:	B