



## Test kit, chlorine (total), model CN-65

**Nr. produs:** 225401  
**RON Preț (fără TVA):** Contact  
**Nu este afișată data expirării**

Produse periculoase



The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specificații

Cei mai mici pași de incrementare:	1
Domeniu de măsurare:	0,2 - 4,0 mg/L & 1 - 20 mg/L
Model test kit:	CN-65
Numar teste:	100
Nume metoda:	Drop count titration/ Thiosulfate
Parametru:	Chlorine, mid and high range - As Total Cl <sub>2</sub>
Platforma:	Drop Count Titration
Smallest increment:	0.2
	1
Tipul carcasei:	D