



SA filter, cone, diameter 110 mm, soil analysis grade, 500/pk

Nr. produs: 2263356

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Designed for sandy soils, SA filters are conically folded and have been optimized for use with the common extraction solutions used in soil nutrient analysis (Mehlich 1, Mehlich 3, Morgan, Modified Morgan, Bray, and Olsen). SA Filters provide the best combination of retention rating, filtration speed and purity requirements for soil nutrient testing (P, K, Ca, Mg). Medium filter retention.

Specificații

Cantitate:	50 /pk
Diametru:	110 mm
Material:	Cellulose