



Test kit, iron color disc, model IR-21

Nr. produs:

2299300

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Iron color disc/TPTZ low range test kit contains 100 tests.

Simple, fast, color disc method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Domeniu de măsurare:	Dual Range
Domeniu de masurare :	0 - 1,2 mg/L Fe
Interval de măsurare standard:	0 - 2.0 mg/L Fe; 0 - 1.2 mg/L Fe
Model:	IR-21
Model test kit:	IR-21
Numar teste:	100
Nume metoda:	Color disc/TPTZ
Parameter:	Iron TPTZ
Parametru:	Iron, low range - As Fe**
Platforma:	Disc culoare
Smallest increment:	0.01
	0.05
Tip:	Dublu metric
Tipul carcasei:	D