



## Test kit, ascorbic acid, model ASC-1

**Nr. produs:**

**2308100**

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

**Produse periculoase**



Model ASC-1 Ascorbic Acid Test Kit was developed in cooperation with beverage producers and uses the drop-count titration method to measure ascorbic acid (vitamin C). Utilizes drop count titration with a range of 10-200 mg/L. Contains approximately 100 tests.

**Utilizes drop count titration - a simple, economical titration method**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Specificații

Domeniu de măsurare:	10 - 200 mg/L
Model test kit:	ASC-1
Nota de subsol:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Numar teste:	100
Nume metoda:	Drop count titration/  Iodometric
Parametru:	Ascorbic Acid
Platforma:	Drop Count Titration
Smallest increment:	10
Tipul carcasei:	D