



A-1 broth medium, dehydrated, 500 g

Nr. produs:

2309934

Produse periculoase

ARTICOL ÎNVECHIT

Acest articol nu mai este disponibil.



A-1 Broth Medium, Dehydrated. Broth for 24-hour MPN method for fecal coliforms. See Microbiological Testing for prepared medium. Bottle/500 g.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Specificații

Cantitate: 500 g

Media Type: A-1 Broth Medium