



Filter, acrocap, pore size 0.2 µm

Nr. produs: 2323810
RON Preț (fără TVA): Contact
Nu este afișată data expedierii

For sterilization or clarification of up to three liters of culture media or other aqueous solutions. Ideal for sterilization of solutions that tend to foam when filtered under vacuum. Inlet accepts 6.4 mm (0.25") tubing or most male luer fittings.

Specificații

Diametru:	NA
Dimensiune por:	0.2 µm
Filter Type:	Filter, acrocap
Packing:	NA