



Filter holder, for vacuum filtration, glass

Nr. produs: 234000

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Beaker-style funnel (300 mL) with 47 mm membrane support for filtration of aqueous samples. Fritted borosilicate glass membrane support for even filtration across surface of membrane. Funnel attaches evenly to base via aluminum clamp. Filter flask sold separately.

Specificații

Cantitate:	1
Conținutul cutiei:	Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8).
Diametru:	47 mm

Conținutul cutiei

Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8).