



Test kit, oxygen scavenger, color disc, model OS-1, 50 tests

Nr. produs: 2349500

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Produse periculoase



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 50 tests.

Simple and fast color disc method

Measures common oxygen scavengers: DEHA, hydroquinone, ISA, and MEKO

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Domeniu de măsurare:	0 - 52 µg/L DEHA (and other Oxygen Scavengers)
Domeniu de masurare :	0 - 1360 µg/L DEHA (and other Oxygen Scavengers)
Model:	OS-1
Model test kit:	OS-1
Numar teste:	50
Nume metoda:	Color disc/Iron reduction
Parametru:	Oxygen Scavenger - As DEHA, hydroquinone, erythorbic acid, and methylethylketoxime
Platforma:	Disc culoare
Tipul carcasei:	R