



Test kit, nitrogen, ammonia, model NI-SA

Nr. produs: 2428700
RON Preț (fără TVA): Contact
 Nu este afișată data expirării

Produse periculoase



Designed for testing ammonia in freshwater or seawater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

Designed for testing ammonia in freshwater or seawater

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specificații

Domeniu de măsurare:	0 - 2.4 mg/L NH ₃ -N
Model test kit:	NI - SA
Numar teste:	100
Nume metoda:	Color disc/Salicylate
Parametru:	Nitrogen, Ammonia, mid range - For freshwater or seawater**
Platforma:	Disc culoare
Smallest increment:	0.1
Tipul carcasei:	D