



Reagent set, molybdenum, LR, 50 mL

Nr. produs:

2440200

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



For determination of Low Range Molybdate Molybdenum by the Ternary Complex method (Hach method 8169). Range: 0.02 - 3.00 mg/L as Mo. Sample Size: 50 mL. Set contains: Molybdenum 1 Reagent Powder Pillows and Molybdenum 2 Reagent Solution for approximately 100 tests.

Specificații

Domeniu de măsurare:	0.02 - 3.00 mg/L Mo
Domeniu de masurare :	50 mL
Instrument:	DR/3000, DR/2000
Numar metoda:	8169
Numar teste:	100
Nume metoda:	Ternary Complex