



Reagent set, hardness (total), buret, HR

Nr. produs:

2447700

Produse periculoase

ARTICOL ÎNVECHIT

Acest articol nu mai este disponibil.



Reagent set for the determination of Total Hardness by the ManVer® 2 Buret Titration Method. Range: 100-25,000 mg/L as CaCO₃. Contains: Buffer Solution (Hardness 1), ManVer® 2 Hardness Indicator Powder Pillows and 0.200 N EDTA (TitraVer®) Titrant for up to approximately 100 tests.

Specificații

Caracteristici speciale:	Reagent set, Buret, HR
Domeniu de măsurare:	100 - 25000 mg/L CaCO ₃
EPA compliant:	Da
Numar metoda:	8222
Numar teste:	100
Nume metoda:	EDTA Titration
Parametru:	Hardness (Total)