



## Hardness (Calcium & Total Sequential) HR Reagent Set, Digital Titrator

**Nr. produs:**

**2448700**

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



For sequential determination of Calcium and Total Hardness by EDTA Digital Titrator titration. Hach method 8329. Range: 100 - 4000 mg/L CaCO<sub>3</sub>. Set contains Buffer Solution (Hardness 1), CalVer 2 Calcium Indicator Powder Pillows, ManVer 2 Hardness Indicator Powder Pillows, Potassium Hydroxide Solution (8 N), Sulfuric Acid (5.25 N) and 0.800 M EDTA Digital Titrator Cartridge. Approximately 40 - 100 tests.

### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### Versatile

More than 50 parameter/range combinations are available

### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

---

## Specificații

|                          |                                         |
|--------------------------|-----------------------------------------|
| Caracteristici speciale: | Reagent set, HR, Digital Titrator       |
| Domeniu de măsurare:     | 100 - 4000 mg/L CaCO <sub>3</sub>       |
| Numar metoda:            | 8329                                    |
| Numar teste:             | 40 - 100                                |
| Nume metoda:             | EDTA Titration                          |
| Parametru:               | Hardness (Calcium and Total Sequential) |
| Platforma:               | Titration digital                       |