



Test kit, phenols, model PL-1

Nr. produs: 2483600

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Produse periculoase



Utilizing the simple color disc method with a range of 0 - 4 mg/L, this kit contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Domeniu de măsurare:	0 - 4.0 mg/L
Interval de măsurare standard:	0 - 4.0 mg/L
Model:	PL-1
Model test kit:	PL-1
Numar teste:	100
Nume metoda:	Color disc/
	4-aminoantipyrine
Parameter:	Phenols
Parametru:	Phenols
Platforma:	Disc culoare
Smallest increment:	0.2 mg/L
Tipul carcasei:	D