



## Flask, Erlenmeyer, Wide Mouth, 2000 mL, Subdivisions 600 mL

**Nr. produs:** 2489454

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Graduated glass flask with wide mouth for easier titrations. Reinforced tops with round finish provide maximum strength. Durable graduations indicate approx. volume. Stopper not included. Requires Size 13 rubber/silicone stopper.

---

### Specificații

|             |                  |
|-------------|------------------|
| Capacitate: | 2000 mL          |
| Material:   | Glass Wide Mouth |
| Unitate:    | 1 each           |