



Filter, membrane, black, pore size 0.45 µm, diameter 47 mm, 200/pk, w/out pads

Nr. produs: 2490300
RON Preț (fără TVA): Contact
Nu este afișată data expirării

Provides contrast for microbiological and particulate analysis. Modified polyethersulfone creates dark background for excellent contrast useful in enumerating opaque colonies and light-colored particulate. White grid line aids enumeration. Certified for membrane filtration analysis. Sterile and individually packaged. Pore size 0.45 µm. Diameter: 47 mm.

Specificații

Cantitate:	200 /pk
Diametru:	47 mm
Dimensiune por:	0.45 µm
Filter Type:	Filter
Material:	Modified polyethersulfone