



Sample Cell: 1" Square Glass 10 mL matched set of 8

Nr. produs:	2495408
RON Preț (fără TVA):	Contact
Disponibil	

Accessories for spectrophotometers/spectral measurements of liquids.

Accessories for spectrophotometers/spectral measurements of liquids.

Comes as a matched set of 8 cells.

Use in Combination with Hach Spectrophotometers

Hach spectrophotometers, chemistries, and methods are designed, developed, and manufactured to optimize your results.

Specificații

Cantitate:	8 /pk
Nota de subsol:	1 pair of cells is included with purchase of DR 5000 UV-Vis, DR 3800, DR 2800, or DR 2700 Spectrophotometer